# WAVE ABSORPTION IN SOLUTIONS OF ORGANIC DYES. ANALYSIS OF SOLUTION COMPOSITION

#### <u>Equipment</u>

- 1. Spectrophotometer
- 2. Solution of fluorescein (FL)  $2 \times 10^{-4} \text{ mol/dm}^3 (200 \,\mu\text{M})$
- 3. Solution of Bengali rose (BR) 2 x  $10^{-4}$  mol/dm<sup>3</sup> (200  $\mu$ M)

#### **Course of the practical**

#### I. PREPARATION OF DYE SOLUTIONS.

- 1. You are given water solutions of two dyes: fluorescein (FL) and Bengali rose (BR). The concentration of both stock solutions is  $2 \times 10^{-4} \text{ mol/dm}^3$  (200  $\mu$ M).
- 2. Dilute stock solution of each dye with distilled water so that you obtain the series of solutions (concentrations 2, 3, 4, 5, 6, 7 and 8  $\mu$ mol/dm<sup>3</sup>). Total volume of each solution is 10 cm<sup>3</sup>. Mix all the solutions carefully. (Put a rubber stopper into the test tube and invert the tube several times. Before transferring the stopper to the next solution dry the stopper with a paper tissue).

#### **II. FINDING MOLAR ABSORPTION COEFFICIENTS OF STUDIED DYES**

#### A. Measurements

Turn on the spectrophotometer with the switch on its back. Wait at least 15 minutes. Use this time to prepare dye solutions.

## ATTENTION 1: keep the spectrophotometer's chamber closed while turning the device on.

## ATTENTION 2: the cuvettes can be touched only at side (opaque) walls.

# ATTENTION 3: do not throw away FL and BR solutions (8 µmol/dm<sup>3</sup>) after the measurements. They are necessary for preparing the dye mixture.

- 1. Place a cuvette filled with distilled water inside the spectrophotometer. Check the position of the cuvette a small triangle on it should overlap a white line inside the measurement chamber. Close the lid.
- 2. Set the measurement mode to Absorption using "A/T/C" button. You should see a letter A in the lower right corner of the display.
- 3. Using the figures of absorption spectra find the analytical wavelengths (with accuracy of 1 nm) for the monomer of fluorescein (FL) and Bengali rose (BR) and write them down to the Table 1.
- 4. Set the wavelength chosen in the point II.3 for FL using the buttons "100", "010" and "001" for setting hundreds, tens and units, respectively.
- 5. Accept the wavelength with the button "BLANK".
- 6. Remove the cuvette with distilled water from the spectrophotometer.
- 7. Fill the second, dry cuvette with the solution of FL of the lowest concentration. Place the cuvette into the spectrophotometer's chamber. Close the lid. After a few seconds read off the absorption intensity from the display and write it down in the Table 2.

- 8. Pour the solution off the cuvette, dry it by pressing against a paper towel. Fill the cuvette with the FL of the next higher concentration and repeat point 7.
- 9. Perform the analogous measurements for Bengali rose.

# **B.** Graphical presentation of the results

- 1. Based on the data from the Table 2 plot the dependence  $A(c) = \varepsilon \cdot l \cdot c$  (absorbance A versus the concentration c), for both dyes.
- 2. Based on the above plots calculate molar absorption coefficients ( $\epsilon$ ) for both dyes. Write the obtained values down to the boxes below the Table 2.

# III. ANALYSIS OF MIXTURE OF DYES.

- In a clean test tube prepare dye mixture by mixing two different volumes of FL and BR solutions (concentration 8 µmol/dm<sup>3</sup>). Use higher volume of FL solution and lower volume of BR solution. Write these volumes down to the Table 3.
- 2. The analysis of absorption spectra of both dyes gives us some important pieces of information:
  - a. absorption of FL solution is equal to 0 in the wavelength range 530 560 nm,
  - b. for the wavelength range between ca. 450 nm and ca. 520 nm the absorption of a mixture is a sum of absorptions of FL and BR solutions.

# Taking into account the additivity of the absorption, we can find the concentration of RB in the dye mixture using the calibration curve prepared for RB in point II.B.1.

- 3. Finding BR concentration in the dye mixture using the calibration plot. Since absorption of FL solution is equal to 0 in the wavelength range 530 560 nm, the concentration of BR in the mixture ( $x_{BR}$ ) can be read off the calibration curve prepared previously for BR (point II.B.1).
  - a. measure the absorbance of the dye mixture at the wavelength of maximal absorption for RB and write it down in Table 4.
  - b. mark this absorbance value on the calibration plot obtained in point II.B.1,
  - c. read the BR concentration in dye mixture off this plot and write it down to the Table 4.
- 4. Knowing the volumes and concentrations of FL and BR solutions used to prepare a mixture, calculate theoretical FL and BR concentrations ( $x_{FL}$ <sup>theoret</sup> and  $x_{BR}$ <sup>theoret</sup>) in the mixture. Write these concentrations down to the Table 5. The value  $x_{BR}$ <sup>theoret</sup> should approximately be the same as the value of  $x_{BR}$  estimated according to the procedure described in the point III.2.

## **Required theoretical knowledge:**

- 1. What is spectroscopic analytical method?
- 2. Types of chemical bonds in organic compounds and names of molecular orbitals.
- 3. What is a chromophore in a structure of organic molecule.
- 4. Ground and excited state of a molecule, mechanism of excitation.
- 5. Types of electronic transitions in molecules.
- 6. What is UV-VIS spectroscopy?
- 7. Light absorption laws:
  - a) Lambert Law (I absorption law)
  - b) Lambert-Beer law (II absorption law)
  - c) absorption addition law (III absorption law)
- 8. What is electronic absorption spectrum and what are its parameters ?

- 9. What is a monomer and what is an aggregate of organic dye in a water solution.
- 10.Condition, at which molecules of organic dye in a solution exist in a monomeric phase.
- 11. Main reasons of deviations from Lambert-Beer law.
- 12. Spectral analysis of mixture of two dyes:
  - a) interpretation of equation about additivity of absorption of two organic dyes
  - b) describe the method of determination of unknown concentration of and bengali rose (RB) in a mixture

Wrocław Medical University Department of Biophysics and Neuroscience	Practical No 36 Wave absorption in solutions of organic dyes. Analysis of solution composition	
		Faculty:
Student names		Date:
Grade:	Tutor's signature:	

# I. Analysis of dyes' absorption spectra

# Table 1. Analytical wavelengths

$\lambda_{max} (nm)$	FL	$\lambda_{max} (nm)$	BR

# II. Graphical presentation of the results

Table 2. Dependence of absorbance on dye concentration

No.	FL concentration mol/dm <sup>3</sup>	A (FL)	BR concentration mol/dm <sup>3</sup>	A (BR)
1.				
2.				
3.				
4.				
5.				
6.				
7.				

$\varepsilon_{BR} =$
$\epsilon_{\sf BR}$ =

# III. Analysis of mixture of dyes

Table 3. Dyes' volumes used to prepare the mixture

FL volume	BR volume

Table 4. Finding BR concentration in the dye mixture using calibration curve

λ (nm)	Α(λ)	X <sub>BR</sub>

Table 5. Theoretical FL and BR concentrations in dye mixture.

$$x_{FL}^{theoret} =$$
  $x_{BR}^{theoret} =$